DCT (SP8) manual

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A. MatLabo

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Command flow:

```
dct_start [gtSetup, gtPreprocessing, gtSegmentationGUI]
```

```
gtSegmentationGUI
```

```
dct_analysis [gtSetupIndexing, gtSetupForwardSimulation, gtSetupReconstruction]
```

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Check carefully that the following are correct:

collection_dir	/ssd1/masayuki/DCT/230418017/L1_230418017	char	Collection directory
name	L1_230418017	char	Name of the dataset
dir	/ssd1/masayuki/DCT/230418017/L1_230418017	char	Directory in which to analyse the data
xdet	2048	double	Detector ROI size X or U (raw image size in pixels)
ydet	2048	double	Detector ROI size Y or V (raw image size in pixels)
nproj	1800	double	Number of images in *180 DEGREES* of scan
refon	3600	double	References after how many images
nref	30	double	How many reference images in a group
ndark	1	double	How many dark images taken
pixelsize	0.0031	double	Detector pixelsize (mm/pixel)
energy	30	double	Beam energy (keV)
dist	8	double	Sample-detector distance (mm) [computed]
sensortype	frelon	char	Camera type ('frelon'/'kodak4mv1'/'marana')
type	360degree	char	DCT scan type ('360degree','180degree', etc)
interlaced_turns	0	double	Interlaced scan? 0 for normal scan, 1 for one extra turn, etc
mono_tune	0	double	Monochromator was tuned after N reference groups, or 0 for not tuned
rotation_axis	vertical	char	Rotation axis orientation ('vertical'/'horizontal')
distortion	none	char	Distortion correction file with path (or 'none')
no_direct_beam	<input type="checkbox"/>	logical	Special scan with no direct beam (taper frelon, offset detector)?
detector_definition	inline	char	Definition of the detector type ('inline'/'vertical')
nof_phases	1	double	Number of phases in sample to be analysed
flip_images	<input type="checkbox"/>	logical	Do you want to flip the images left-right for some reason?
online	<input type="checkbox"/>	logical	Is the analysis online?
scan_type	ftomo	char	Acquisition macro: {'ebs_tomo','fscan_v1','finterlaced'}

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Figure 2: gtModifyStructure Version 003@pool2

Edit the values, click "OK" when you are happy

Parameters dealing with arbitrary geometry:

Parameter	Value	Type	Description
beamdir	1 0 0	double	Beam direction in LAB reference (unit row vector)
rotpos	0 0 0	double	Rotation axis position (arbitrary point on axis) in LAB
deflabX	Along the beam direction.	char	Description how Lab X direction was chosen [for records only]
deflabY	Right-handed from Y=cross(Z,X).	char	Description how Lab Y direction was chosen [for records only]
deflabZ	Along rotation axis. Positive away from sample stage.	char	Description how Lab Z direction was chosen [for records only]
labunit	mm	char	LAB units (default is mm) [for records only]
rotdir	0 0 1	double	Rotation axis direction in LAB (unit row vector); omega is right-handed rotation

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Number of crystallographic phases to be analysed in this sample? : [1] >>

Figure 1: gtModifyStructure Version 003@pool2

Edit the values, click "OK" when you are happy

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
Do you want to load it from a .cif file? [y/n] : [n] >>

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Do you want to load reflection data from .csv or .dat files? (recommended) [y/n] : [n] >>

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Do you want to look at the contents of cryst for this phase? [y/n] : [y] >>

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Start moving and correction process now? [y/n] : [y] >>

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Do you want to verify/change the (default) preprocessing parameters? [y/n] : [n] >>

Do you want to redefine the direct beam BoundingBox (currently [583 738 770 910])? [y/n] : [n] >>

Do you want to (re)-calculate the Rotation Axis position (currently 1024.5)? [y/n] : [n] >>

Do you want to (re)-define the Sample BoundingBox (currently [655 697 740 774])? [y/n] : [n] >>

Do you want to (re)-define the Active Area of the detector? [y/n] : [n] >>

Do you want to perform sample drift analysis for this measurement? [y/n] : [n] >>

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Do you want to reset all the parameters for segmentation? [y/n] : [n] >>

Recompile functions for OAR? Not normally needed! [y/n] : [n] >>

Check intensities in direct beam ? [y/n] : [n] >>

Drop and recreate fullmedianvals table? [y/n] : [y] >>

Offset the median value of each full image to 0 (recommended)? [y/n] : [y] >>

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Drop and recreate difblob tables? [y/n] : [y] >>

Drop and recreate difspot tables? [y/n] : [y] >>

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

(上)solve_conflict_with_intensity = false (下) = true

GtGrainsAnalysis : L1_DCT@pool2

Limits: [-1, 609] (456, 269, 223) Auto min Zoom Flip L/R Flip U/D XY YZ ZX
Phase: 1, Grain 0 clim max YX ZY XZ Info

Slice: 223

50
100
150
200
250
Y

6

ID
 Rand
 Volu...
 Co...
 Caxis
 IPF
 Rvec
 Sch...
 Mos...
 Ann...

Set Cmap Options
Supported Cmap
Open Cmap Editor
Save Cmap Figure
Update Cmap

Calculate SchmidF
Calculate Mosaicity
Pole Figure
Inverse Pole Figure

Rand:
0.000

Select/Deselect
Reconstruct & Seg...
Deselected Grains

Phase 01
Complete Volume
Dilated Volume
Assemble & D.
Auto Assemble

Phase Volume
Complete Volume
Dilated Volume

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Options for IPF colormap:

sampleDir	<input type="text" value="0 0 1"/>	double	Sample direction in the sample CS for IPF calculation
phaseid	<input type="text" value="1"/>	double	Current phase
saturate	<input checked="" type="checkbox"/>	logical	Saturate colors in the IPF map
save	<input type="checkbox"/>	logical	Save colormap on disk?
background	<input checked="" type="checkbox"/>	logical	Do we want the background color in the cmap?

Clear Quit Check **OK**

Slide through the slices

Status: Done. Phase ID (0 = All): 1 Apply Grain ID (0 = All): 0 Apply Open